



A.B.N. 88 007 386 156  
**MEDICAL  
SUPPLIES PTY. LTD.**

Factory 8, 21 Jersey Road  
(P.O. Box 184)  
Bayswater, Vic, 3153

Phone: 03 .9720 6177  
Fax: 03 .9720 6178

Email: info@rgmedical.com.au  
www.rgmedical.com.au

## **ICS PACIFIC G135L MEDICAL** Medical Refrigerator

### **ICS PACIFIC G135L Medical Under counter or bench top Vaccine Fridge**

### **Fully compliant with Quality Care Pharmacy Program (QCPP) Cold Chain Standards**

#### **Features**

- Internal Illumination
- High/low temperature display with reset button
- Electronic temperature control and display
- Audible temperature alarm system
- External alarm facility
- Maintenance free condenser
- Ventilated cooling system
- Automatic defrost
- Door lock
- Stackable
- Self-closing door
- 4 Adjustable shelves plus base
- Capacity: 155 Litre
- 3 year Manufacturers Warranty

#### **Temperature Range:**

+2°C to + 8°C @ 30°C ambient temp  
60% humidity

#### **Dimensions:**

**Width:**  
600mm

**Depth:**  
610mm

**Height:**  
850mm



**Brand New stock  
(not scratch and dent)**

**Sale price \$1199 +GST**

**Scratch & Dent Clearance**

**Sale Price \$980 +GST**

Save 40% off RRP (\$1650 +GST)

All fridges come with 36 month Manufacturers Warranty  
1x Foam marks on door & stain on bulkhead  
1x Dent on top LHS corner  
1x Minor scratch on top  
1x Small Dent on RHS

Freight charges will apply